




FLT3 Antibody

Product Code	CSB-PA008343
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P36888
Immunogen	Synthesized peptide derived from Human Flt3 around the non-phosphorylation site of Y599.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	FLT3; CD135; FLK2; STK1; Receptor-type tyrosine-protein kinase FLT3; FL cytokine receptor; Fetal liver kinase-2; FLK-2; Fms-like tyrosine kinase 3; FLT-3; Stem cell tyrosine kinase 1; STK-1; CD antigen CD135
Immunogen Species	Homo sapiens (Human)
Target Names	FLT3

Image



Western Blot analysis of HEPG2 cells using Flt3 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.