



# ARHGEF9 Antibody

|                            |  |
|----------------------------|--|
| <b>Product Code</b>        | CSB-PA008344   |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.  |
| <b>Uniprot No.</b>         | O43307   |
| <b>Immunogen</b>           | Synthesized peptide derived from the C-terminal region of Human ARHGEF9.   |
| <b>Raised In</b>           | Rabbit   |
| <b>Species Reactivity</b>  | Human, Mouse, Rat  |
| <b>Tested Applications</b> | IHC, ELISA; IHC:1:100-1:300, ELISA:1:40000   |
| <b>Form</b>                | Liquid   |
| <b>Conjugate</b>           | Non-conjugated   |
| <b>Storage Buffer</b>      | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.  |
| <b>Purification Method</b> | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.                            |
| <b>Isotype</b>             | IgG  |
| <b>Alias</b>               | ARHGEF9; ARHGEF9; KIAA0424; Rho guanine nucleotide exchange factor 9; Collybistin; PEM-2 homolog; Rac/Cdc42 guanine nucleotide exchange factor 9 |
| <b>Immunogen Species</b>   | Homo sapiens (Human)   |
| <b>Target Names</b>        | ARHGEF9  |