




ADNP Antibody

Product Code	CSB-PA008363
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H2P0
Immunogen	Synthesized peptide derived from the N-terminal region of Human ADNP.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ADNP; ADNP1; KIAA0784; Activity-dependent neuroprotector homeobox protein; Activity-dependent neuroprotective protein
Immunogen Species	Homo sapiens (Human)
Target Names	ADNP

Image



Western Blot analysis of LOVO cells using ADNP Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.