



HEY2 Antibody

Product Code	CSB-PA008379
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UBP5
Immunogen	Synthesized peptide derived from the Internal region of Human HRT2.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	HEY2; BHLHB32; CHF1; GRL; HERP; HERP1; HRT2; Hairy/enhancer-of-split related with YRPW motif protein 2; Cardiovascular helix-loop-helix factor 1; hCHF1; Class B basic helix-loop-helix protein 32; bHLHb32; HES-related repressor protein 2; Ha
Immunogen Species	Homo sapiens (Human)
Target Names	HEY2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.