



SMUG1 Antibody

Product Code	CSB-PA008405
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q53HV7
Immunogen	Synthesized peptide derived from the Internal region of Human SMUG1.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat, Monkey
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	SMUG1; Single-strand selective monofunctional uracil DNA glycosylase
Immunogen Species	Homo sapiens (Human)
Target Names	SMUG1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.