



PLD3 Antibody

Product Code	CSB-PA008426
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IV08
Immunogen	Synthesized peptide derived from the Internal region of Human PC-PLD3.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	IHC, ELISA; IHC:1:100-1:300, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	PLD3; Phospholipase D3; PLD 3; Choline phosphatase 3; HindIII K4L homolog; Hu-K4; Phosphatidylcholine-hydrolyzing phospholipase D3
Immunogen Species	Homo sapiens (Human)
Target Names	PLD3