




POTEH/POTEG Antibody

Product Code	CSB-PA008442
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6S545/Q6S5H5
Immunogen	Synthesized peptide derived from the Internal region of Human POTE-14/22.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	POTEH; A26C3; ACTBL1; POTE22; POTE ankyrin domain family member H; ANKRD26-like family C member 3; Prostate; ovary, testis-expressed protein on chromosome 22; POTE-22; POTEG; A26C2; POTE14; POTE ankyrin domain family member G; ANKRD26-like
Immunogen Species	Homo sapiens (Human)
Target Names	POTEH/POTEG

Image



Western Blot analysis of AD293 using POTE-14/22 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.