



Phospho-ALK (Y1604) Antibody

Product Code	CSB-PA008446
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UM73
Immunogen	Synthesized peptide derived from Human ALK around the phosphorylation site of Y1604.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ALK; ALK tyrosine kinase receptor; Anaplastic lymphoma kinase; CD antigen CD246
Immunogen Species	Homo sapiens (Human)
Target Names	ALK
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.