



Phospho-ALOX5 (S272) Antibody

Product Code	CSB-PA008472
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P09917
Immunogen	Synthesized peptide derived from Human 5-LO around the phosphorylation site of S272.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ALOX5; LOG5; Arachidonate 5-lipoxygenase; 5-LO; 5-lipoxygenase
Immunogen Species	Homo sapiens (Human)
Target Names	ALOX5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.