



PRPF6 Antibody

Product Code	CSB-PA008476
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O94906
Immunogen	Synthesized peptide derived from the C-terminal region of Human PRPF6.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	IHC, ELISA; IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	PRPF6; C20orf14; Pre-mRNA-processing factor 6; Androgen receptor N-terminal domain-transactivating protein 1; ANT-1; PRPF6 homolog; U5 snRNP-associated 102 kDa protein; U5-102 kDa protein
Immunogen Species	Homo sapiens (Human)
Target Names	PRPF6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.