



GAST Antibody

Product Code	CSB-PA008504
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01350
Immunogen	Synthesized peptide derived from the Internal region of Human Gastrin.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	GAST; GAS; Gastrin
Immunogen Species	Homo sapiens (Human)
Target Names	GAST
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.