



FUT2 Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA008506 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q10981 |
| Immunogen | Synthesized peptide derived from the Internal region of Human FUT2. |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | IHC, ELISA; IHC:1:100-1:300, ELISA:1:20000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Purification Method | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Isotype | IgG |
| Alias | FUT2; SEC2; Galactoside 2-alpha-L-fucosyltransferase 2; Alpha(1; 2)FT 2; Fucosyltransferase 2; GDP-L-fucose:beta-D-galactoside 2-alpha-L-fucosyltransferase 2; SE2; Secretor blood group alpha-2-fucosyltransferase; Secretor factor; Se |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | FUT2 |