





Phospho-FYN (Y530) Antibody

Product Code	CSB-PA008519
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P06241
Immunogen	Synthesized peptide derived from Human Fyn around the phosphorylation site of Y530.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	FYN; Tyrosine-protein kinase Fyn; Proto-oncogene Syn; Proto-oncogene c-Fyn; Src-like kinase; SLK; p59-Fyn
Immunogen Species	Homo sapiens (Human)
Target Names	FYN


Image



Western Blot analysis of LOVO+UV cells using Phospho-Fyn (Y530) Polyclonal Antibody



Western blot analysis of 3T3 lysis using Phospho-Fyn (Y530) antibody.



Western blot analysis of 3T3 lysis using Phospho-Fyn (Y530) antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.