



FZD2 Antibody

Product Code	CSB-PA008522
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14332
Immunogen	Synthesized peptide derived from the Internal region of Human Frizzled-2.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	FZD2; Frizzled-2; Fz-2; hFz2; FzE2
Immunogen Species	Homo sapiens (Human)
Target Names	FZD2