



PLEKHG4 Antibody

Product Code	CSB-PA008584
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q58EX7
Immunogen	Synthesized peptide derived from the Internal region of Human Puratrophin 1.
Raised In	Rabbit
Species Reactivity	Human,Monkey
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	PLEKHG4; PRTPHN1; Puratrophin-1; Pleckstrin homology domain-containing family G member 4; PH domain-containing family G member 4; Purkinje cell atrophy-associated protein 1
Immunogen Species	Homo sapiens (Human)
Target Names	PLEKHG4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.