



# DAK Antibody

|                     |   |
|---------------------|---|
| Product Code        | CSB-PA008649  |
| Storage             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.   |
| Uniprot No.         | Q3LXA3  |
| Immunogen           | Synthesized peptide derived from the Internal region of Human DHA Kinase.   |
| Raised In           | Rabbit  |
| Species Reactivity  | Human,Mouse,Rat   |
| Tested Applications | WB, ELISA; WB:1:500-1:2000, ELISA:1:20000   |
| Form                | Liquid  |
| Conjugate           | Non-conjugated  |
| Storage Buffer      | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.   |
| Purification Method | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Isotype             | IgG   |
| Alias               | DAK; Bifunctional ATP-dependent dihydroxyacetone kinase/FAD-AMP lyase; cyclizing                                      |
| Immunogen Species   | Homo sapiens (Human)  |
| Target Names        | TKFC  |