



STEAP2 Antibody

Product Code	CSB-PA008693
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NFT2
Immunogen	Synthesized peptide derived from the C-terminal region of Human STEAP2.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	IHC, IF, ELISA; IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	STEAP2; PCANAP1; STAMP1; Metalloreductase STEAP2; Prostate cancer-associated protein 1; Protein up-regulated in metastatic prostate cancer; PUMPCn; Six-transmembrane epithelial antigen of prostate 2; SixTransMembrane protein of prostate 1
Immunogen Species	Homo sapiens (Human)
Target Names	STEAP2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.