



Phospho-GATA3 (S308) Antibody

Product Code	CSB-PA008717
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P23771
Immunogen	Synthesized peptide derived from Human GATA-3 around the phosphorylation site of S308.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	IF, ELISA; IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	GATA3; Trans-acting T-cell-specific transcription factor GATA-3; GATA-binding factor 3
Immunogen Species	Homo sapiens (Human)
Target Names	GATA3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.