




GATA3/GATA2 Antibody

Product Code	CSB-PA008718
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P23771/P23769
Immunogen	Synthesized peptide derived from Human GATA-3 around the non-phosphorylation site of S308.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IHC, IF, ELISA; IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	GATA3; Trans-acting T-cell-specific transcription factor GATA-3; GATA-binding factor 3; GATA2; Endothelial transcription factor GATA-2; GATA-binding protein 2
Immunogen Species	Homo sapiens (Human)

Image



Western blot analysis of Mouse-brain lysis using GATA-2/3 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.