



GCHFR Antibody

Product Code	CSB-PA008753
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P30047
Immunogen	Synthesized peptide derived from the Internal region of Human GFRP.
Raised In	Rabbit
Species Reactivity	Human,Rat
Tested Applications	IHC, ELISA; IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	GCHFR; GFRP; GTP cyclohydrolase 1 feedback regulatory protein; GFRP; GTP cyclohydrolase I feedback regulatory protein; p35
Immunogen Species	Homo sapiens (Human)
Target Names	GCHFR
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.