



OR8B2/OR8B3 Antibody

Product Code	CSB-PA008775
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96RD0/Q8NNG8
Immunogen	Synthesized peptide derived from the Internal region of Human Olfactory receptor 8B2/3.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	OR8B2; Olfactory receptor 8B2; Olfactory receptor OR11-309; OR8B3; Olfactory receptor 8B3; Olfactory receptor OR11-311
Immunogen Species	Homo sapiens (Human)
Target Names	OR8B2/OR8B3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.