



NPAS4 Antibody

Product Code	CSB-PA008779
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IUM7
Immunogen	Synthesized peptide derived from the C-terminal region of Human NPAS4.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	IHC, ELISA; IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	NPAS4; BHLHE79; NXF; PASD10; Neuronal PAS domain-containing protein 4; Neuronal PAS4; Class E basic helix-loop-helix protein 79; bHLHe79; HLH-PAS transcription factor NXF; PAS domain-containing protein 10
Immunogen Species	Homo sapiens (Human)
Target Names	NPAS4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.