



Phospho-GFAP (S38) Antibody

Product Code	CSB-PA008785
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P14136
Immunogen	Synthesized peptide derived from Human GFAP around the phosphorylation site of S38.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	GFAP; Glial fibrillary acidic protein; GFAP
Immunogen Species	Homo sapiens (Human)
Target Names	GFAP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.