




Phospho-HBP1 (S402) Antibody

Product Code	CSB-PA008810
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60381
Immunogen	Synthesized peptide derived from Human HBP1 around the phosphorylation site of S402.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	HBP1; HMG box-containing protein 1; HMG box transcription factor 1; High mobility group box transcription factor 1
Immunogen Species	Homo sapiens (Human)
Target Names	HBP1

Image



Western Blot analysis of A549 cells using Phospho-HBP1 (S402) Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.