



Phospho-GJA1 (S261) Antibody

Product Code	CSB-PA008819
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17302
Immunogen	Synthesized peptide derived from Human Connexin 43 around the phosphorylation site of S261.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	IHC, ELISA; IHC:1:100-1:300, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	GJA1; GJAL; Gap junction alpha-1 protein; Connexin-43; Cx43; Gap junction 43 kDa heart protein
Immunogen Species	Homo sapiens (Human)
Target Names	GJA1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.