




Phospho-PDCD4 (S457) Antibody

Product Code	CSB-PA008870
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q53EL6
Immunogen	Synthesized peptide derived from Human Pdcd-4 around the phosphorylation site of S457.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	PDCD4; H731; Programmed cell death protein 4; Neoplastic transformation inhibitor protein; Nuclear antigen H731-like; Protein 197/15a
Immunogen Species	Homo sapiens (Human)
Target Names	PDCD4

Image



Western Blot analysis of 293 cells using Phospho-Pdcd-4 (S457) Polyclonal Antibody