



GNG5 Antibody

Product Code	CSB-PA008889
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P63218
Immunogen	Synthesized peptide derived from the N-terminal region of Human Gγ 5.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	IF, ELISA; IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	GNG5; GNGT5; Guanine nucleotide-binding protein G(I)/G(S)/G(O) subunit gamma-5
Immunogen Species	Homo sapiens (Human)
Target Names	GNG5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.