



GPR17 Antibody

Product Code	CSB-PA008913
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13304
Immunogen	Synthesized peptide derived from the Internal region of Human GPR17.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	GPR17; Uracil nucleotide/cysteinyll leukotriene receptor; UDP/CysLT receptor; G-protein coupled receptor 17; P2Y-like receptor; R12
Immunogen Species	Homo sapiens (Human)
Target Names	GPR17
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.