




GPS2 Antibody

Product Code	CSB-PA008949
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13227
Immunogen	Synthesized peptide derived from the N-terminal region of Human GPS2.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	GPS2; G protein pathway suppressor 2; GPS-2
Immunogen Species	Homo sapiens (Human)
Target Names	GPS2

Image



Western Blot analysis of RAW264.7 cells using GPS2 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.