




GRIN1 Antibody

Product Code	CSB-PA008974
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q05586
Immunogen	Synthesized peptide derived from Human NMDA ζ 1 around the non-phosphorylation site of S897.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	GRIN1; NMDAR1; Glutamate [NMDA] receptor subunit zeta-1; N-methyl-D-aspartate receptor subunit NR1; NMD-R1
Immunogen Species	Homo sapiens (Human)
Target Names	GRIN1

Image



Western Blot analysis of 3T3 cells using NMDA ζ 1 Polyclonal Antibody