





Phospho-NR3C1 (S203) Antibody

Product Code	CSB-PA008993
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04150
Immunogen	Synthesized peptide derived from Human GR around the phosphorylation site of S203.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	NR3C1; GRL; Glucocorticoid receptor; GR; Nuclear receptor subfamily 3 group C member 1
Immunogen Species	Homo sapiens (Human)
Target Names	NR3C1

Image



Western Blot analysis of A549 cells using Phospho-GR (S203) Polyclonal Antibody



Western Blot analysis of HELA 293T 453 cells using Phospho-GR (S203) Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.