



CXCL3 Antibody

Product Code	CSB-PA009037
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19876
Immunogen	Synthesized peptide derived from the C-terminal region of Human GRO γ .
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	IHC, ELISA; IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CXCL3; GRO3; GROG; SCYB3; C-X-C motif chemokine 3; GRO-gamma(1-73); Growth-regulated protein gamma; GRO-gamma; Macrophage inflammatory protein 2-beta; MIP2-beta
Immunogen Species	Homo sapiens (Human)
Target Names	CXCL3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.