




GSPT1 Antibody

Product Code	CSB-PA009039
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P15170
Immunogen	Synthesized peptide derived from the Internal region of Human eRF3a.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	GSPT1; ERF3A; Eukaryotic peptide chain release factor GTP-binding subunit ERF3A; Eukaryotic peptide chain release factor subunit 3a; eRF3a; G1 to S phase transition protein 1 homolog
Immunogen Species	Homo sapiens (Human)
Target Names	GSPT1

Image



Western Blot analysis of 3T3 cells using eRF3a Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.