




GTF2H2 Antibody

Product Code	CSB-PA009043
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13888
Immunogen	Synthesized peptide derived from the N-terminal region of Human TFIIH p44.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	GTF2H2; BTF2P44; General transcription factor IIH subunit 2; Basic transcription factor 2 44 kDa subunit; BTF2 p44; General transcription factor IIH polypeptide 2; TFIIH basal transcription factor complex p44 subunit
Immunogen Species	Homo sapiens (Human)
Target Names	GTF2H2

Image



Western Blot analysis of COLO205 cells using TFIIH p44 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.