




Phospho-BLNK (Y96) Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA009047 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q8WV28 |
| Immunogen | Synthesized peptide derived from Human BLNK around the phosphorylation site of Y96. |
| Raised In | Rabbit |
| Species Reactivity | Human, Mouse, Monkey |
| Tested Applications | WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:20000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Purification Method | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Isotype | IgG |
| Alias | BLNK; BASH; SLP65; B-cell linker protein; B-cell adapter containing a SH2 domain protein; B-cell adapter containing a Src homology 2 domain protein; Cytoplasmic adapter protein; Src homology 2 domain-containing leukocyte protein of 65 kDa; |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | BLNK |

Image



Western Blot analysis of COLO205 cells using Phospho-BLNK (Y96) Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.