



GPR132 Antibody

Product Code	CSB-PA009059
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UNW8
Immunogen	Synthesized peptide derived from the C-terminal region of Human GPR132.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	GPR132; G2A; Probable G-protein coupled receptor 132; G2 accumulation protein
Immunogen Species	Homo sapiens (Human)
Target Names	GPR132
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.