



Acetyl-H2AFZ (K5) Antibody

Product Code	CSB-PA009100
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0C0S5
Immunogen	Synthesized peptide derived from Human Histone H2A around the acetylation site of K5.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	IHC, ELISA; IHC:1:100-1:300, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	H2AFZ; H2AZ; Histone H2A.Z; H2A/z
Immunogen Species	Homo sapiens (Human)
Target Names	H2AFZ
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.