



Phospho-H3F3A (S31) Antibody

Product Code	CSB-PA009122
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P84243
Immunogen	Synthesized peptide derived from Human Histone H3.3 around the phosphorylation site of S31.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	IHC, IF, ELISA; IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	H3F3A; H3.3A; H3F3; PP781; H3F3B; H3.3B; Histone H3.3
Immunogen Species	Homo sapiens (Human)
Target Names	H3F3A