




Phospho-ACACA (S80) Antibody

Product Code	CSB-PA009195
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13085
Immunogen	Synthesized peptide derived from Human ACC α around the phosphorylation site of S80.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ACACA; ACAC; ACC1; ACCA; Acetyl-CoA carboxylase 1; ACC1; ACC-alpha
Immunogen Species	Homo sapiens (Human)
Target Names	ACACA

Image



Western blot analysis of 3T3 using p-ACC α (S80) antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.