



HLA-DPA1 Antibody

Product Code	CSB-PA009209
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P20036
Immunogen	Synthesized peptide derived from the Internal region of Human HLA-DP α 1.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	HLA-DPA1; HLA-DP1A; HLAB*07:01; HLA class II histocompatibility antigen; DP alpha 1 chain; DP(W3); DP(W4); HLA-SB alpha chain; MHC class II DP3-alpha; MHC class II DPA1
Immunogen Species	Homo sapiens (Human)
Target Names	HLA-DPA1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.