



VN1R5 Antibody

Product Code	CSB-PA009237
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7Z5H4
Immunogen	Synthesized peptide derived from the Internal region of Human V1RL5.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	IF, ELISA; IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	VN1R5; V1RL5; Vomeronasal type-1 receptor 5; G-protein coupled receptor GPCR26; hGPCR26; V1r-like receptor 5
Immunogen Species	Homo sapiens (Human)
Target Names	VN1R5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.