




HOXA5 Antibody

Product Code	CSB-PA009244
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P20719
Immunogen	Synthesized peptide derived from the N-terminal region of Human HoxA5.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	HOXA5; HOX1C; Homeobox protein Hox-A5; Homeobox protein Hox-1C
Immunogen Species	Homo sapiens (Human)
Target Names	HOXA5

Image



Western Blot analysis of 293T cells using HoxA5 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.