



HOXB1/HOXA1/HOXD1 Antibody

Product Code	CSB-PA009246
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P14653/P49639/Q9GZZ0
Immunogen	Synthesized peptide derived from the C-terminal region of Human Hox-A1/B1/D1.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	IHC, ELISA; IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	HOXB1; HOX2I; Homeobox protein Hox-B1; Homeobox protein Hox-2I; HOXA1; HOX1F; Homeobox protein Hox-A1; Homeobox protein Hox-1F; HOXD1; HOX4; HOX4G; Homeobox protein Hox-D1; Homeobox protein Hox-GG
Immunogen Species	Homo sapiens (Human)
Target Names	HOXB1/HOXA1/HOXD1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.