



HOXB2 Antibody

Product Code	CSB-PA009247
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P14652
Immunogen	Synthesized peptide derived from the Internal region of Human HoxB2.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	HOXB2; HOX2H; Homeobox protein Hox-B2; Homeobox protein Hox-2.8; Homeobox protein Hox-2H; K8
Immunogen Species	Homo sapiens (Human)
Target Names	HOXB2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.