




Phospho-HSP90AB1 (S226) Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA009327 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P08238 |
| Immunogen | Synthesized peptide derived from Human HSP90 β around the phosphorylation site of S226. |
| Raised In | Rabbit |
| Species Reactivity | Human, Mouse, Rat |
| Tested Applications | WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:10000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Purification Method | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Isotype | IgG |
| Alias | HSP90AB1; HSP90B; HSPC2; HSPCB; Heat shock protein HSP 90-beta; HSP 90; Heat shock 84 kDa; HSP 84; HSP84 |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | HSP90AB1 |

Image



Western blot analysis of VEC using p-HSP90 β (S226) antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.