



HTR2A Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA009355 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P28223 |
| Immunogen | Synthesized peptide derived from the C-terminal region of Human SR-2A. |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | IHC, IF, ELISA; IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:5000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Purification Method | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Isotype | IgG |
| Alias | HTR2A; HTR2; 5-hydroxytryptamine receptor 2A; 5-HT-2; 5-HT-2A; Serotonin receptor 2A |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | HTR2A |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |