




HTR4 Antibody

Product Code	CSB-PA009360
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13639
Immunogen	Synthesized peptide derived from the Internal region of Human SR-4.
Raised In	Rabbit
Species Reactivity	Human,Rat
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	HTR4; 5-hydroxytryptamine receptor 4; 5-HT-4; 5-HT4; Serotonin receptor 4
Immunogen Species	Homo sapiens (Human)
Target Names	HTR4

Image



Western Blot analysis of 3T3 cells using SR-4 Polyclonal Antibody