



HTR7 Antibody

Product Code	CSB-PA009368
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P34969
Immunogen	Synthesized peptide derived from the C-terminal region of Human SR-7.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	IF, ELISA; IF:1:200-1:1000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	HTR7; 5-hydroxytryptamine receptor 7; 5-HT-7; 5-HT7; 5-HT-X; Serotonin receptor 7
Immunogen Species	Homo sapiens (Human)
Target Names	HTR7