




ICAM3 Antibody

Product Code	CSB-PA009391
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P32942
Immunogen	Synthesized peptide derived from Human CD50 around the non-phosphorylation site of S518.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ICAM3; Intercellular adhesion molecule 3; ICAM-3; CDw50; ICAM-R; CD antigen CD50
Immunogen Species	Homo sapiens (Human)
Target Names	ICAM3

Image



Western Blot analysis of NIH-3T3 cells using CD50 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.