



AKR1CL1 Antibody

Product Code	CSB-PA009427
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5T2L2
Immunogen	Synthesized peptide derived from the Internal region of Human AKR1CL1.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	IF, ELISA; IF:1:200-1:1000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	AKR1CL1; Aldo-keto reductase family 1 member C-like protein 1
Immunogen Species	Homo sapiens (Human)
Target Names	AKR1C8P